

G6PDH LYPHO

(Quantitative/Kinetic Method)

Intended Use

Reagent Kit for quantitative estimation of G6PDH in erythrocytes.

Reagent Provided

1. G6PDH (Co-Enzyme – Substrate)
2. G6PDH (Buffer Reagent)
3. G6PDH (Lysing Reagent)


Diagnostic Significance

G-6-PD deficient subjects are prone to hemolytic episodes and hemolytic anaemia if treated with certain antimalarial or sulfa drugs. The screening for the deficiency is one of the first preventive step to avoid these hemolytic episodes and anaemia. The estimation will indicate deficiency better in males than in females.

Principle

The enzyme Glucose-6-Phosphate Dehydrogenase present in the Red Blood Cells is extracted by lysing the cells using a natural detergent. The extracted enzyme oxidizes Glucose-6-Phosphate to 6-Phosphogluconate and simultaneously reduces co-enzyme NADP to NADPH giving increase in absorbance at 340 nm.

Enzymatic determination of G6PDH activity is based on the following reaction.



Components & Concentration of Reagents

Reagent	Component	Concentration
Active Ingredients	Glucose-6-Phosphate	0.5 mmol/L
	NADPNa ₂	0.75 mmol/L
	Buffer	100 mmol/L
	Detergent	1 mmol/L
	pH7.5±0.1 at 25°C	

Reagent Storage & Stability

G6PDH reagents are stable until the expiry date stated on the label. Reconstituted reagent is stable for 7 days at 2-8°C.

Preparation of Working Reagent

Add 1.1 ml buffer (2G6PDH) to the Vial labeled 1G6PDH. Mix well and use after five minutes.

Specimen Collection & Preservative

Fresh whole blood is the specimen required. Collection of blood by using any one of the anticoagulants such as Citrate, Oxalate or Heparin is recommended.

Preparation of Red Cell Hemolysate

Wash 0.1 ml of whole blood with 2 ml aliquots of physiological saline (0.9%) 3 times and suspend the washed, packed and centrifuged erythrocytes in precooled 0.5 ml of 3G6PDH (lysing

Reagent). Mix well and keep in the refrigerator (2-4°C) for atleast 15 minutes and maximum for 2 hours. Centrifuge the lysate at 3000 r.p.m for 5 minutes prior to use.

Precaution

G6PDH is for IN-VITRO diagnostic use only.

Reagent contains Sodium Azide.

Do Not Ingest

The Lysing Reagent (3G6PDH) should be used cold (0-8°C) to avoid the Decrease in enzyme activity.

Assay Guidelines for Analyzers

Reaction Type	Kinetic
Wavelength	340 nm
Slope of Reaction	Increasing
FlowCell Temperature Setting	30°C
Zero Setting With	Distilled Water
Delay Time	180 sec
Interval Time	60 sec
No. of Readings	4
Working Reagent Volume	1.0 ml
Sample Volume	25 µl (0.025 ml)
Factor	3601 (When Hemoglobin is Measured) 36013 (When RBC Count is measured)
Linearity	Up to (22.5 U/g Hb) or (628 U/10 ¹² RBC'S)

Test Procedure

Determine the Hemoglobin content of the whole blood and the RBC count prior to lysis of the cells.

Pipette Into Test Tubes	Test
Working Reagent	1.0 ml
Hemolysate	0.025 ml

Mix and aspirate. After the initial delay of 180 sec., record the absorbance of test at the interval of one minute for the next 3 minutes at 340 nm. Determine the mean change in absorbance per minute and calculate test results.

Calculation

(i) PDH Activity (U/10¹² RBC)

$$\Delta A/\text{Min} \times \frac{224 \times 10^{12}}{6.22 \times N \times 10^6 \times 1000}$$

$$\Delta A/\text{Min} \times \frac{36013}{N}$$

Where

224 = Total assay volume to sample volume

10¹² = Factor for expressing activity in 10¹² cells

6.22 = Milimolar absorptivity of NADPH at 340 nm.
 $N \times 10^6$ = Number of erythrocytes/cmm
 1000 = Conversion of Cell Count from Count per cmm to count per ml

(ii) G6PDH Activity (U/gHb)

$$\Delta A/\text{Min} \times \frac{224 \times 100}{6.22 \times \text{Hb}(\text{g/dl})}$$

$$\Delta A/\text{Min} \times \frac{3601}{\text{Hb}(\text{g/dl})}$$

Where

100 = Factor to convert to 100 ml
 224 = Total assay volume to sample volume
 6.22 = Milimolar absorptivity of NADPH at 340 nm.
 Hb(g/dl) = Hemoglobin Concentration

Example

Assay of specimen which had a red cell count of $4.6 \times 10^6/\text{mm}^3$ and a hemoglobin concentration of 15.2 g/dl result in $\Delta A/\text{min}$ at 30°C of 0.028.

$$\text{G6PDH}(\text{U}/10^{12} \text{RBC}) = \frac{0.028 \times 36013}{4.6} = 219$$

$$\text{G6PDH}(\text{U/gHb}) = \frac{0.028 \times 3601}{15.2} = 6.63$$

Normal Values

G6PDH Activity (at 30°C) 146-376 U/ 10^{12} RBCs
 4.6 – 13.5 U/g Hb

Linearity

The maximum G6PDH activity which may be measured by this procedure is approx: (22.5 U/gHb) or (628 U/ 10^{12} RBC's). If the value is greater than the reported value repeat determination using 5 μl blood as sample and multiply results by 2.

Notes

If the G6PDH activity is very low, measure the absorbance change for 5 minutes after addition of buffered substrate and divide by 5 to obtain $\Delta A/\text{minute}$ and calculate the test results. Use the formula given under "TEST RESULTS".

References

1. Kachmar J.F. Moss S.W. Enzymes. In Fundamentals of Clinical Chemistry Ed. By N.W. 2. Teitz. Saunders Philadelphia. 1976 p.p 666-672.
2. Gross R.T, Hurwitz R.E., Marks P.A. Red Cell Glucose-6-Phosphat dehydrogenase deficiency. J. Clin. Invest (1978) 37.176.

Symbols

IVD

In Vitro Diagnostics



Caution

LOT

Batch No.

CONT

Content



Read Instructions



Storage Temperature

REF

Catalogue No.



Product Expiry Date



Manufactured By



Date of Manufacture



Keep Dry



Fragile



Keep away from sun light



KEE DIAGNOSTICS PVT LTD

CIN: U24231DL2004PTC128343

A-8, 3RD FLOOR, NARAINA INDUSTRIAL AREA, PHASE - II, NEW DELHI 110028

W: www.keediagnostics.com

E: info@keediagnostics.com

T: 011 43136000



ISO 9001: 2015 & ISO 13485: 2016 Certified Company