

LIPASE

(Colorimetric Method)


Intended Use

Lipase is a reagent set for determination of Lipase activity in serum and plasma based on colorimetric method.

Principle

The chromogenic lipase substrate 1,2-O-dilauryl-rac-glycero-3-glutaric acid-(6-methylresorufin) ester is cleaved by the catalytic action of alkaline lipase solution to form 1,2-O-dilauryl-rac-glycerol and an unstable intermediate glutaric acid-(6-methylresorufin)ester. This decomposes spontaneously in alkaline medium to form glutaric acid and methylresorufin.

The lipase activity in the specimen is proportional to the formation of methylresorufin in the reaction and can be determined photometrically at 578 nm.



Reagents provided

- 1. Reagent R₁
- 2. Reagent R₂

Reagent storage and stability

The reagent kit should be stored at 2 - 8°C and is stable till the expiry date indicated on the label. R₁ and R₂ reagents are stable till expiry at 2 - 8°C. Do not freeze the reagent.

Specimen collection and preservation

Blood should be collected in a clean dry container. Do not use calcium complexing anticoagulants such as EDTA, citrate and fluoride. Lipase in serum/plasma is stable for 1 week at 15 - 25°C/1 week at 2 - 8°C/1 year at - 20°C.. Centrifuge samples containing precipitate before performing the assay. The samples should be brought to room temperature prior to use.

Assay guidelines for Analyzer

Reaction type	Fixed Time
Reaction slope	Increasing
Wavelength	578 nm
Flow cell temperature	37°C
Zero setting	Distilled Water
Delay Time	60 seconds
Read Time	60 seconds
Reagent Volume	1000µl (800µl + 200 µl)
Sample volume	10 µl (0.010 ml)
Linearity	250 U/L

Perform the assay as given below:

	Calibrator	Sample
Reagent 1	800 µl	800 µl
Reagent 2	200µl (0.200 ml)	200 µl (0.200 ml)
Sample	10 µl (0.010 ml)	10 µl (0.010 ml)

1. Mix thoroughly and transfer the assay mixture immediately to the thermo stated cuvette and start the stop watch simultaneously.
2. Record the first reading at 60th second and subsequently 1 more readings with 60 seconds interval at 578 nm.

Calculation

$$\text{Activity of Lipase in U/L} = \frac{\Delta \text{Abs. /min of Sample}}{\Delta \text{Abs. /min of Calibrator}} \times \text{Conc. of Calibrator}$$

Normal Range

Adults: upto 60 U/L (0.22-1.00 μ kat/L)

Note: Expected range varies from population to population. It is therefore recommended that each laboratory should establish its own normal range.

Limitations

- If the Lipase activity exceeds 250 U/L, dilute the specimen with normal saline and repeat the assay. In such cases the results obtained should be multiplied with the dilution factor to obtain correct Lipase activity.

Quality Control

To ensure adequate quality control, it is recommended that each batch should include normal and abnormal commercial reference control serum. It should be realized that the use of quality control material checks both instrument and reagent functions together. Factors which might affect the performance of this test include proper instrument function, temperature control, cleanliness of glassware, and accuracy of pipetting.

References

1. Greiling H, Gressner AM, eds, Lehrbuch der Klinischen **Chemie und Pathobiochemie**, 3rd ed. Stuttgart / New York: Schattauer Verlag 1995.
2. Keller H, ed. **Klinisch-chemische Labordiagnostik für die Praxis**, 2nd ed. Stuttgart/ New York: Georg Thieme Verlag 1991:354-361.
3. Tietz NW, Shuey DF. Lipase in serum - the elusive enzyme: An overview. **Clin Chem** 1993;39 (5): 746-756.
4. **Tietz NW, ed. Clinical Guide to Laboratory Tests**, 3rd ed. Philadelphia PA: WB Saunders Company 1995:865.
5. **Neumann U, Junius M, Batz HG, et al.** New substrates for the optical determination of lipase; EP 2072521987.

Symbols

IVD	In Vitro Diagnostics.	 Caution.	 Keep away from sun light.	 Date of Manufacture.
LOT	Batch No.	 Read Instructions.	 Fragile.	 Product Expiry Date.
CONT	Content.	 Storage Temperature.	 Keep Dry.	 Manufactured By.
REF	Catalogue No.			



KEE DIAGNOSTICS PVT LTD

CIN: U24231DL2004PTC128343 .

A subsidiary of KEE PHARMA LTD

A-8 . 3RD FLOOR. NARAINA INDUSTRIAL AREA, PHASE 2.
NEW DELHI-110028.

W. www.keediagnostics.com

E. info@keediagnostics.com, T. 011-43136000.



ISO 9001: 2015 & ISO 13485:2016 Certified Company