

MICROALBUMIN

(Turbilatex Method)

Intended Use

Microalbumin is a reagent kit used for the determination of albumin (urine) concentration in urine. The reagents are for *in-vitro* diagnostic use.

Principle

The kit utilizes latex-enhanced immunoturbidimetry to measure the mALB level in human serum or plasma. During the test, mALB in the sample binds with the specific mALB antibody which is coated on latex particles to cause agglutination. The turbidity caused by agglutination is detected optically by chemistry analyzer. The change in absorbance is proportional to the level of mALB in the sample. The actual concentration is obtained by comparing with a calibration curve with known concentrations.

Reagent provided

1. Reagent 1 – Buffer
2. Reagent 2 – Latex
3. Reagent 3 – Calibrator

Working Reagent Preparation

Prepare working reagent by mixing **Reagent R1** and **Reagent R2** in the ratio 4:1 as per the number of tests required.

Reagent storage and stability

Store all the reagents at 2-8°C.

Specimen collection and preservation

Follow standard laboratory procedures to collect urine samples and store them at 2-4°C for up to 2 days or at -20°C for up to 1 months. Avoid repeated freezing and thawing.

Serial Dilution Step

	1st	2nd	3rd	4th	5th
Calibrator	100 µl	50 µl from 1st Tube	50 µl from 2nd Tube	50 µl from 3rd Tube	50 µl from 4th Tube
Normal Saline	0	50 µl	50 µl	50 µl	50 µl
Ratio of Dilution	Neat	1/2	1/4	1/8	1/16

System Parameters

Method	Fixed Time (2-Point)
Wavelength	630 nm
Zero Setting	Distilled Water
Temperature Setting	37°C
Incubation Temperature	37°C
Incubation Time	-
Delay Time	10 sec
Read Time	120 sec
No. of Reading	2
Reagent 1 Volume	800 µl
Reagent 2 Volume	200 µl
Sample Volume	10 µl
Standard Concentration	Refer Calibrator Vial
Units	mg/L
Factor	-
Reaction Slope	Increasing
Linearity	300 mg/L

Assay guidelines for Analyzers

These reagents may be used both for manual assay and in several automatic analyzers.

Reagent	Calibrator (C)	Test (T)
R1 Buffer	800 µl	800 µl
R2 Antibody	200 µl	200 µl
Bring up the temperature of determination. Then add		
Calibrator	10 µl	-
Sample	-	10 µl

Mix well, after about 10 sec (37°C) read the absorbance A1 of the test (T) and calibrator (C) against **Distilled water**. After exactly 120 sec (for all temperature) read the absorbance A2 of the test (T) and calibrator (C). Calculate $\Delta A/\text{min}$ (A2-A1) for the test and calibrator.

Calculation

$$\text{mALB Conc. (mg/L)} = \frac{\text{Abs. Test}}{\text{Abs Cal}} \times \text{Conc. Cal}$$

Normal Range

Up to 25 mg/dL or 30 mg/day

It is recommended that each laboratory establish its own normal range.

Limitation

1. For in vitro diagnostic use only, Do not pipette by mouth. Avoid contact with skin and mucous membranes.
2. All the calibrators and controls must be considered as human & animal sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.
3. Exercise the normal precautions required for handling laboratory reagents.
4. Reagents with different lot numbers should not be interchanged or mixed.

Quality Control

To ensure adequate quality control, it is recommended that each batch should include a commercially available control material with established values determined by this method. It should be realized that the use of quality control material checks both instrument and reagent functions together. Factors which might affect the performance of this test include proper instrument function, temperature control, and cleanliness of glassware and accuracy of pipetting.

References

1. Persijn, J.P., et al, Clin. Acta 35:91, (1971).
2. Stookey, L.L., Anal. Chem. 42:779, (1970).
3. Tietz, N.W., Fundamentals of Clinical Chemistry Philadelphia, W.B. 1. Saunders, pp. 923-929, (1976).
4. Weissman, N., Pileggi, V.J., in Clinical Chemistry: Principles and Technics,
5. 2nd Ed., R.J. Henry et al. editors, Hagerstown (MD), Harper & Row, pp. 692-693, (1974).
6. Young, D.S. et al, Clin. Chem. 21:1D, (1975).
7. Henry, J.B., Clinical Diagnosis and Management by Laboratory Methods,
8. Philadelphia, W.B. Saunders, p. 1434, (1984).

Symbols

IVD	In Vitro Diagnostics.		Caution.		Keep away from sun light.		Date of Manufacture.
LOT	Batch No.		Read Instructions.		Fragile.		Product Expiry Date.
CONT	Content.		Storage Temperature.		Keep Dry.		Manufactured By.
REF	Catalogue No.						



KEE DIAGNOSTICS PVT LTD

CIN: U24231DL2004PTC128343

A subsidiary of KEE PHARMA LTD

A-8 . 3RD FLOOR. NARAINA INDUSTRIAL AREA, PHASE 2.

NEW DELHI-110028.

W. www.keediagnosics.com.

E. info@keediagnosics.com, T. 011 43136000

CE ISO 9001 : 2015 & ISO 13485 : 2016 Certified Company